

SPECIFICATIONS

Components	REF	W634P0024	W634P0028	W634P0025	W634P0027
Sealed Pouch(pcs)		1	2	5	20
Extraction Buffer		1	2	5	20
Disposable Sterile Swab(pcs)		1	2	5	20
Waste Bag(pcs)		1	2	5	20
Instruction for Use(pcs)		1	1	1	1

STORAGE AND STABILITY

- The test kit should be stored at 2–30°C (storage in refrigerator is permitted). Do not store the kit in the freezer. Improper storage may result in an inaccurate result.
- The test cassette is sensitive to humidity and temperature. Once removed from foil pouch, test cassette is stable for up to 1 hour.
- The test kit is stable until the expiration date printed on outer package. Do not use it beyond the expiration date.
- The test cassette must remain in the sealed pouch until use.

LIMITATIONS OF PROCEDURE

- This reagent is designed to detect 2019-nCoV antigen in human nasal swab specimen.
- Failure to follow the instructions for the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
- Reading the test results earlier than 15 minutes or later than 20 minutes may give incorrect results.
- The sample collection process will affect the accuracy of the test, such as improper sample collection, improper sample storage, etc.
- This reagent is a qualitative assay. As it is with any diagnostic procedure, a confirmed 2019-nCoV infection diagnosis should only be made by a physician after evaluating all clinical and laboratory findings.
- Negative test results may occur if the level of antigen in a sample is below the detection limit of the test, or from improper sample collection, and the negative results are not intended to exclude other non 2019-nCoV virus infections.
- A negative test result does not rule out a coronavirus infection and does not exempt you from the applicable rules for spread control (e.g. contact restrictions and protective measures). Symptomatic patient must seek further testing (e.g. immediate PCR), even though negative result occurs.

Wondfo Quick Guide

2019-nCoV Antigen Test (Lateral Flow Method)



Operation Video

WHAT DOES THE KIT TEST?


Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is a rapid test that is used for laymen of detecting novel coronaviruses (2019-nCoV) N protein antigen extracted from the nasal swab specimen. It is intended as an aid in the diagnosis of coronavirus infection disease (COVID-19) for asymptomatic patients and/or symptomatic patients within 7 days after onset of symptoms, which is caused by 2019-nCoV.

For *in vitro* diagnostic use only. For self-testing use.

According to usability study on laymen user, the test can be correctly performed for anyone age over 18. However, children under the age of 18 should be supported or under the supervision by an adult.

MAKE SURE YOUR TEST KIT CONTAINS

- Sealed Pouch
- Extraction Buffer (400μL)
- Disposable Sterile Swab
- Waste Bag
- Instruction for Use
- Tube Rack (in the outer box)



Other required items (not included in the test kit)

Clock or timer

- Positive test results do not exclude co-infections with other pathogens or identify specific 2019-nCoV virus subtype (like SARS-CoV virus), and cannot necessarily determine whether a person is infectious.
- SARS-CoV virus variants including Delta (B.1.617.2) and Omicron (B.1.1.529) have been detected out by Wondfo 2019-nCoV Antigen Test (Lateral Flow Method).

QUESTION & ANSWER

Q1. How does the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) work?

The Wondfo 2019-nCoV Antigen Test is an antigen test that is to detect the presence of protein fragments (antigen) from the 2019-nCoV in nasal swab specimen.

Q2. What is the difference between a COVID-19 antigen, molecular, and antibody test?

There are different kinds of tests for diagnosing COVID-19. Molecular tests (also known as PCR tests) detect genetic material from the virus. The Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is an antigen test which detects small parts or proteins from the virus. Antigen tests are very specific for the virus but are not as sensitive as molecular tests. Another type of test is an antibody test. A COVID-19 antibody test detects antibodies that have been made by your immune system in response to a previous COVID-19 infection.

Q3. Will this test hurt?

No, the disposable sterile swab is not sharp and it should not hurt. Sometimes the swab can feel slightly uncomfortable or tickly.

Q4. Why do I swab both nostrils?

Swabbing both nostrils gives you the best chance of collecting sufficient sample to generate an accurate result. It has been observed in some cases that only one nostril has detectable virus, so it is important to collect from both nostrils. Correct swabbing is important to obtain a correct result.

Q5. What does it mean if I have a positive test result?

A positive result means that you may have COVID-19 disease. Please contact your doctor for further medical suggestion. It is likely you will be asked to isolate yourself at home to avoid spreading the virus to others, wear a face mask when recommended and wash your hands regularly with soap and water. A positive result does not in any way guarantee that you are or will be immune and therefore cannot (or can no longer) become infected.

Q6. What does it mean if I have a negative test result?

A negative result means the virus that causes COVID-19 was not found in your sample. A negative test result does not guarantee that you do not or have never

had COVID-19, nor does it confirm whether or not you are currently contagious.

Do you have cold symptoms in addition to the negative at-home test? Since the at-home test does not provide complete certainty, you should assume that you have COVID-19. You can contact your doctor to find out if another test is needed. In the meantime, try to avoid leaving your home and have as little contact as possible with others, including the people you live with. Use disposable tissues and throw them straight in the bin. Sneeze and cough into the crook of your elbow. Wash your hands regularly and wear a face mask. Are your symptoms getting worse (difficulty breathing, high fever, etc.)? Contact your doctor/health provider immediately.

Q7. How accurate is the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)?

The test has been shown in field clinical evaluations performed by professional health care persons to correctly identify 99.90% (963 out of 964) of 2019-nCoV negative samples (known as the test's specificity), when compared with RT-PCR test. Further, in field clinical evaluations conducted in Germany, US and UK, the test correctly identified 98.53% (335/340) 2019-nCoV negative samples when performed by lay users, and compared with RT-PCR test. The test has also been shown in field clinical evaluations performed by professional health care persons to correctly identify 92.40% (304 out of 329) of 2019-nCoV positive samples (known as the test's sensitivity), when compared with RT-PCR test. Further, in field clinical evaluations conducted in Germany, US and UK, the test correctly identified 92.96% (66/71) of 2019-nCoV positive samples when performed by lay users, and compared with RT-PCR test.

Among the 319 RT-PCR positive samples, when the CT≤25, the detectable rate of the test is 100% (145 out of 145); when the CT≤28, the detectable rate of the test is 98.02% (198 out of 202); when the CT≤30, the detectable rate of the test is 95.88% (233 out of 243); when the CT≤35, the detectable rate of the test is 92.83% (285 out of 307); Thereinto 468 are from symptomatic patients. The test can correctly identify 100% (161 out of 161) 2019-nCoV negative samples and 91.67% (242 out of 264) 2019-nCoV positive samples, when the specimens collected within 7 days after onset of symptoms.

Thereinto 425 are from asymptomatic patients. The test can correctly identify 100% (380 out of 380) 2019-nCoV negative samples and 95.56% (43 out of 45) 2019-nCoV positive samples, for the specimens collected from asymptomatic patients.

Detection limit

The detection limit of Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is 5.0×10³ TCID₅₀/mL.

Q8. Is there any chance that I get a "false" negative result with this test?

It is possible for this test to give an incorrect negative (false negative) result". This means that you could still have COVID-19 even though the test result is negative. If your result is negative and you still experience

symptoms related to COVID-19, such as fever, cough and/or shortness of breath, you should seek help from your healthcare provider.

Q9. Is there any chance that I get an incorrect positive result?

There is a very small chance that this test gives you a positive result that is incorrect (false positive). If you get a positive result, you should self-isolate and seek medical help from you healthcare provider.

Q10. I have used the test but no colored band appears at control line (C). What should I do?

If there is no colored band appears at control line (C) within 15 minutes of performing the test, then the test has not worked. You should test again, using a new test, taking care to follow the instruction.

Q11. Can any medication or medical conditions affect the results?

Yes, it may affect your test result, consult your doctor, and always read the medication manufacturers' instructions for any medication you are taking before conducting the test.

Besides, the test result will not be affected by the presence of following potentially interfering and cross-reactive substances in the specimens:

Mucin	Human coronavirus 229E	Enterovirus
Chloraseptic (Menthol/Benzocaine)	Human coronavirus OC43	Respiratory syncytial virus
Naso GEL (NeilMed)	Human coronavirus NL63	Rhinovirus Type 1A
CVS Nasal Drops (Phenylephrine)	MERS-coronavirus	Haemophilus influenzae Type b
Afrin (Oxymetazoline)	Human Adenovirus 1	Streptococcus pneumonia
CVS Nasal Spray (Cromolyn)	Human Metapneumovirus 3 (hMPV-3) Type B1	Streptococcus pyogenes
Zicam	Parainfluenza virus Type 1	Candida albicans
Homeopathic (Alkalol)	Parainfluenza virus Type 2	pooled human nasal wash
Sore Throat Phenol Spray	Parainfluenza virus Type 3	Bordetella pertussis
Tobramycin	Parainfluenza virus Type 4A	Mycoplasma pneumonia
Mupirocin	Influenza A (H3N2)	Chlamydia pneumonia
Fluticasone Propionate	Influenza A (H1N1)	Legionella pneumophila
Tamiflu (Osetamivir Phosphate)	Influenza B (Victoria lineage)	Staphylococcus aureus
	Influenza B (Yamagata lineage)	Staphylococcus epidermidis

Q12. What are the possible risks of this test?

Possible Risks:


- Discomfort during the sampling
- Incorrect test results (see Interpreting Results and Limitations Sections).

BIBLIOGRAPHY

- Centers for Disease Control and Prevention (CDC). Interim Guidelines for Collecting, Handling, and Testing for Patients with Suspected Novel Influenza A (H1N1) Virus Infection. Available online at: <https://www.cdc.gov/h1n1flu/specimencollection.htm>
- Song F, Zhang X, Zha Y, Liu W. COVID-19: Recommended sampling sites


Step 8

Unscrew the small cap at the top of the Extraction Buffer Tube.




Step 9

Lay the Cassette flat and add 4 drops processed specimen into the sample well. Wait for 15 minutes and read the results. **DO NOT** read results after 20 minutes.



Step 10

The used device may be disposed of with normal household waste in accordance with the applicable local regulations.



Step 11

Re-apply hand sanitizer.



Step 12

Result interpretation

Positive (+)

Two red lines will appear. One on the top half and one on the bottom half. COVID-19 was detected.

NOTE: It does not matter if one of the lines that make up the test line (T) is lighter or darker than the other; the result is "Positive". (Please refer to Q5 in the Instruction for Use for details)

at different stages of the disease. J Med Virol. 2020;92(9):1383-1385. doi:10.1002/jmv.25892.

- Tu YP, O'Leary TJ. Testing for Severe Acute Respiratory Syndrome-Coronavirus 2: Challenges in Getting Good Specimens, Choosing the Right Test, and Interpreting the Results. Crit Care Med. 2020;48(11):1680-1689. doi:10.1097/CCM.0000000000004594.

INDEX OF SYMBOL

	Do Not Reuse		See Instruction for Use		Expiry Date
	Manufacturing Date		Keep Dry		Batch Number
	Keep Away from Sunlight		Manufacturer		Catalog #
	Tests Per Kit		In Vitro Diagnostic Use		Authorized Representative
	Store Between 2~30°C				

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City, Luogang District, 510663 Guangzhou, P.R.China
Tel: 0086-20-3229-9890/ 0086-20-3229-9786
Website: www.wondfo.com.cn
E-mail: service@wondfo.com.cn

CE **0123**
Qarad EC-REP BV
Pas 257
2440 Geel
Belgium

Suppliers of disposable sterile swab


- Miracleclean Technology Co., Ltd.** **0197**(according to Directive 93/42/EEC)
No. 18, Rongshuxia Industrial Zone Tongle Community, Longgang District Shenzhen 518116 Guangdong China
EC representative name: Share Info Consultant Service LLC Repräsentanzbüro
EC representative address: Heerdter Lohweg 83, 40549 Düsseldorf, Germany
- Jiangsu Changfeng Medical Industry Co., Ltd.** **0197** (according to Directive 93/42/EEC)
Touqiao Town, Guangling District Yangzhou 225109 Jiangsu P.R.China
EC representative name: Lins Service & Consulting GmbH
EC representative address: Obere Seegasse 34/2, 69124 Heidelberg, Germany
- Medico Biomedical Technology Co., Ltd.** **0413** (according to Directive 93/42/EEC)
Room 201 of Building 14th and Building 17th, Hengyi Lane, Yuanhu Road, Zhangbei Industrial Park, Longcheng Street, Longgang District, Shenzhen, Guangdong, China
EC representative name: Welikang Ltd
EC representative address: Enterprise Hub, NW Business Complex, 1 BeraghmoreRd., Derry, BT488SE, N. Ireland, UK
- Jiangsu Hanheng Medical Technology co., Ltd.** **0197** (according to Directive 93/42/EEC)
16-B4, #1 North Qingyang Road, Tianning District, 213017 Changzhou, Jiangsu, China
EC representative name: Luxus Lebenswelt GmbH
EC representative address: Kochstr. 1, 47877, Willich, Germany
- Shenzhen KangDaAn Biological Technology Co., Ltd.** **0197** (according to Directive 93/42/EEC)
East-1, 3rd floor, Building 2, Shunheda Factory Liuxiandong industrial zone, Xili street Nanshan district, Shenzhen 518055 Guangdong P.R. China
EC representative name: Share Info Consultant Service LLC Repräsentanzbüro
EC representative address: Heerdter Lohweg 83, 40549 Düsseldorf, Germany

Negative Result (-)

A single red line on the top half. COVID-19 was not detected. (Please refer to Q6 in the Instruction for Use for details)

Invalid

If you see no line appearing on the top half, the test is invalid. It is recommended to repeat the test from collecting a new nasal swab.



Rev. A4 Rel.:2022/02/08