


Wondfo®



+ NEW POCT SOLUTION



PRODUCT CATALOG | 2020



WONDFO PROFILE **Stock Code: 300482**

Guangzhou Wondfo Biotech Co., Ltd., located in Guangzhou, China, is a globally recognized fast growing POCT company devoted to the R&D, manufacture and sales of rapid diagnostic reagents and related equipment in IVD (*in vitro* diagnostic) field.

Founded in 1992, the company is continually expanding its scale. Apart from the Guangzhou headquarters, Wondfo also owns two overseas R&D center, one branch in Chicago and more than ten overseas offices. Wondfo went public in June, 2015 (Stock Code: 300482), and has more than 1,500 employees.

As one of the earliest high-tech biological companies focusing on rapid diagnostics for toxicology, infectious diseases, fertility, cardiovascular disease, tumor, respiratory disease, renal injury and diabetes. Wondfo has won good reputation around the world. In addition, Wondfo has got CE certificates, FDA certificates, FSC certificates, ISO 9001:2015 and ISO 13485:2016 quality system certification.

We are dedicated to making the world a better place by providing people with high quality products and services.



01 Professional R&D Team

Wondfo has been consistently paying attention to R&D and regarding innovation as the core of its development. As a result, there are highly educated employees in the company. And it owns 4 well-equipped research laboratories all over the world.

02 Strong Production Capacity

With the international leading production lines and highly experienced staff, its daily production capacity reaches to 1.5 million test strips. Services are also tailored to professional distributors and partnering affiliates.

03 Integrated Services System

Covering more than 180 countries and areas, Wondfo provides an exciting opportunity to distribute quality products at very competitive price. Every aspect of the company's operation is strictly measured, which provides customers with all-around technical support and counselling all the time.

Focus on *POCT* Industry

Cardiovascular




Respiratory



Diabetes



Renal Injury



8

Major fields

Fertility



Tumor



DOA



Infection



Applications



Pharmacy
Drugstore



Laboratory



Hospital
Clinic
GP Office




Anti-drug

CERTIFICATES



ISO 13485



ISO9001



CE



FDA



Guangzhou Wondfo Biotech Co., Ltd.

CONTENTS

Qualitative Rapid Diagnostic Family

Fertility Tests 02

hCG Pregnancy Test
Sperm Concentration Test
LH Ovulation Test
FSH Urine Test

Infectious Disease Tests 03

SARS-CoV-2 Ag/Ab Test
HBsAg Hepatitis B Surface Antigen Test
HBV Hepatitis B Virus Test
HBsAb Hepatitis B Surface Antibody Test
HCV Hepatitis C Virus Test
HIV 1/2 Blood Test (2 Lines)
HIV 1/2 Blood Test (3 Lines)
HIV 1/2 Oral Mucosal Transudate Test
Malaria P.f Test
Malaria P.f/Pan Test
Malaria P.f/P.v Test
H.Pylori Test
Flu A Test
Flu B Test
Flu A&B Test
Dengue NS1 Antigen Test
Dengue IgG/IgM Antibody Test
Combined Dengue NS1 Antigen and Dengue IgG/IgM Antibody Test
HAV Hepatitis A Virus Test
HEV Hepatitis E Virus Test
HIV-HCV-HBsAg-Syphilis Test
HIV-HCV-HBsAg-HBcAb-Syphilis Test
Syphilis Test
Syphilis-HIV Test
Chlamydia Test
Adenovirus Test
Rotavirus Test
Strep A Test

Tumor Marker Tests 04

FOB Fecal Occult Blood Test
PSA Prostate Specific Antigen Test

Drug of Abuse Tests (Urine) 04-06

AMP Amphetamine Test
BAR Barbiturates Test
CAF Caffeine Test
BZO Benzodiazepines Test
COC Cocaine Test
MET Methamphetamine Test
MDMA Test
BUP Buprenorphine Test
MOP Morphine Test
MTD Methadone Test
OPI Opiate Test
PCP Phencyclidine Test
TCA Tricyclic Antidepressants Test
THC Marijuana Test
OXY Oxycodone Test
KET Ketamine Test
PPX Propoxyphene Test
EDDP Test
TRA Tramadol Test
K2 Synthesis Cannabis Test
FTY Fentanyl Test
COT Cotinine Test
ETG Ethylglucuronide Test
MQL Methaqualone Test

Multi-Drug Tests 06

Multi-Drug 6 Test
Multi-Drug T-Cube
Multi-Test Panel
Multi-Test T-Cup
Multi-Test Q-Cup
Multi-Test K-Cup
Multi-Test T-Cup Plus
Multi-Test Q-Cup Plus

Other Rapid Tests 06

Urinalysis Reagents Strip
Troponin I Test
Alcohol Saliva Test

Quantitative Detection Family 07-11

Dry Chemistry Analyzer 07

Blood Gas Analyzer 08

Optical Coagulation Analyzer 09

Blood Glucose Monitoring System 09

Quantitative Detection Family 10

Immunofluorescence Quantitative Test System 11

Qualitative Rapid Diagnostic Family

- Infectious Disease Tests
- Tumor Marker Tests
- Alcohol Tests
- Cardiac Tests
- Urinalysis Reagents Strip



• Fertility Tests



Pregnancy & Ovulation tests



▲ DM-102 (Disposable)

Digital Pregnancy Test



Sperm Concentration Test Kit

• DOA Urine Tests



Cassette

Strip



Panel



Q-Cup

T-Cup

Key-Cup

• DOA Urine Tests



T-Cup Plus (Applicable with Vacuum Tube)

• DOA Oral Fluid Tests



DOA Oral Fluid Tests

Fertility Test: Pregnancy Tests

Cat. No.	Test Item	Format	Width	Specimen	Sensitivity	Qualification
W1-S	hCG Pregnancy Test	Strip	2.5/3.0/3.5/4.0 mm	Urine	25/20/10 mIU/mL	CE CLIA WAIVED 0123
W1-C	hCG Pregnancy Test	Cassette	4.0 mm	Urine	25/20/10 mIU/mL	CE CLIA WAIVED 0123
W1-C I	hCG Pregnancy Test	Cassette	3.5 mm	Urine	25/20/10 mIU/mL	CE CLIA WAIVED 0123
W1-M	hCG Pregnancy Test	Midstream	6.0 mm	Urine	25/20/10 mIU/mL	CE CLIA WAIVED 0123
W1-M II	hCG Pregnancy Test	Midstream	5.0 mm	Urine	25/20/10 mIU/mL	CE CLIA WAIVED 0123
W1-M III	hCG Pregnancy Test	Midstream	6.0 mm	Urine	25/20/10 mIU/mL	CE CLIA WAIVED 0123
W1-M4	hCG Pregnancy Test	Midstream	3.5 mm	Urine	25/20/10 mIU/mL	CE CLIA WAIVED 0123
W1-M6	hCG Pregnancy Test	Midstream	6.0 mm	Urine	25/20/10 mIU/mL	CE CLIA WAIVED 0123
W1-S3P	hCG Pregnancy Test	Strip	2.5/3.0 mm	Urine/Serum	25 mIU/mL	CE CLIA WAIVED
W1-C3P	hCG Pregnancy Test	Cassette	4.0 mm	Urine/Serum	25 mIU/mL	CE CLIA WAIVED
DM-102	Digital Pregnancy Test (Disposable)	Midstream	6.0 mm	Urine	25 mIU/mL	CE CLIA WAIVED 0123
W121-S	Sperm Concentration Test	Strip	2.5/3.0/3.5 mm	Semen	15x10 ⁶ /mL	CE 0123
W121-C	Sperm Concentration Test	Cassette	4.0 mm	Semen	15x10 ⁶ /mL	CE 0123

Fertility Tests: Ovulation Tests

Cat. No.	Test Item	Format	Width	Specimen	Sensitivity	Qualification
W2-S	LH Ovulation Test	Strip	2.5/3.0/3.5 mm	Urine	25 mIU/mL	CE CLIA WAIVED 0123
W2-C	LH Ovulation Test	Cassette	4.0 mm	Urine	25 mIU/mL	CE CLIA WAIVED 0123
W2-C I	LH Ovulation Test	Cassette	3.5 mm	Urine	25 mIU/mL	CE CLIA WAIVED 0123
W2-M	LH Ovulation Test	Midstream	6.0 mm	Urine	25 mIU/mL	CE CLIA WAIVED 0123
W2-M II	LH Ovulation Test	Midstream	5.0 mm	Urine	25 mIU/mL	CE CLIA WAIVED 0123
W2-M4	LH Ovulation Test	Midstream	3.5 mm	Urine	25 mIU/mL	CE CLIA WAIVED 0123
W2-M6	LH Ovulation Test	Midstream	6.0 mm	Urine	25 mIU/mL	CE CLIA WAIVED 0123
W52-S	FSH Urine Test	Strip	2.5/3.0/3.5 mm	Urine	25 mIU/mL	CE CLIA WAIVED 0123
W52-C	FSH Urine Test	Cassette	4.0 mm	Urine	25 mIU/mL	CE 0123
W52-M	FSH Urine Test	Midstream	6.0 mm	Urine	25 mIU/mL	CE 0123

Infectious Disease Tests

Cat. No.	Category	Test Item	Specimen	Format	Sensitivity	Qualification
W195	COVID-19 <small>NEW</small>	SARS-CoV-2 Ab	S/P/WB	Cassette	/	CE CFDA
W003-C		HBsAg	S/P/WB	Cassette	4 IU/mL	
W003-S		HBsAg	S/P/WB	Strip	4 IU/mL	
W004-C		HBsAb	S/P/WB	Cassette	30 mIU/mL	
W004-S		HBsAb	S/P/WB	Strip	30 mIU/mL	
W040-P		HBV(HBsAg, HBsAb, HBeAg, HBeAb, HBcAb)	S/P/WB	Panel	/	
W005-C	Hepatitis	HCV	S/P/WB	Cassette	2 NCU/mL	
W005-S		HCV	S/P/WB	Strip	2 NCU/mL	
W071-C		HAV-IgM	S/P/WB	Cassette	/	CE
W072-C		HEV-IgM	S/P/WB	Cassette	/	CE
W026-C		HIV-HCV-HBsAg-Syphilis Test	S/P/WB	Cassette	/	
W61-P		HIV-HCV-HBsAg-HBcAb-Syphilis Test	S/P/WB	Cassette	/	
W006-C4P2		HIV 1/2(2 Lines)	S/P/WB	Cassette	/	PQ
W752-C4G3	HIV	HIV 1/2(2 Lines)	S/P/WB	Cassette	/	CE <small>0123</small>
W050-C		HIV 1/2(3 Lines)	S/P/WB	Cassette	/	
W063-C		HIV 1/2 Oral Mucosal Transudate Test	Oral Mucosal Transudate	Cassette	/	
W107		Syphilis-HIV	Whole Blood	Cassette	2 NCU/mL	CE
W34-S4P		Syphilis	S/P/WB	Strip	2 NCU/mL	CE
W34-C4P		Syphilis	S/P/WB	Cassette	2 NCU/mL	CE
W35-C	STD	Chlamydia	Swab/Urine	Cassette	5×10 ⁴ IFU/mL	CE <small>0123</small>
W036		Gonorrhea	Swab	Cassette	1×10 ⁵ CFU/mL	CE
W37-C		Malaria HRP2(P.f)	Whole Blood	Cassette	100 p/μL	CE
W37-S		Malaria HRP2(P.f)	Whole Blood	Strip	100 p/μL	CE
W056-C	Malaria	Malaria HRP2/pLDH(P.f/P.v)	Whole Blood	Cassette	100 p/μL	CE
W62-C		Malaria HRP2/pLDH(P.f/Pan)	Whole Blood	Cassette	100 p/μL	CE
W028-C		H.Pylori Antigen	Feces	Cassette	/	CE
W38-C		H.Pylori Antibody	Serum/Plasma	Cassette	/	CE
W38-S	H.Pylori	H.Pylori Antibody	Serum/Plasma	Strip	/	CE
W71-C		Flu A&B	Nasal Swab/Throat Swab	Cassette	Flu A: 5.1×10 ⁵ TCID50/mL Flu B: 1.5×10 ⁶ TCID50/mL	CE <small>0123</small>
W72-C		Flu A	Nasal Swab/Throat Swab	Cassette	Flu A: 5.1×10 ⁵ TCID50/mL	CE <small>0123</small>
W73-C	Flu	Flu B	Nasal Swab/Throat Swab	Cassette	Flu B: 1.5×10 ⁶ TCID50/mL	CE <small>0123</small>
W109-C		Dengue IgG/IgM Antibody Test	S/P/WB	Cassette	/	CE
W041-C		Dengue NS1 Antigen Test	S/P/WB	Cassette	/	CE
W111	Dengue	Combined Dengue NS1 Antigen and IgG/IgM Antibody Test	S/P/WB	Cassette	/	CE
W112-C		Rotavirus Test	Feces	Cassette	/	CE
W133-C		Adenovirus Test	Feces	Cassette	/	CE
W39-C	Others	Strep A	Swab	Cassette	/	CE <small>0123</small> CLIA WAIVED

Tumor Marker Tests

Cat. No.	Test Item	Format	Width	Specimen	Sensitivity	Qualification
W31-S	PSA	Strip	2.5/3.0 mm	Serum/Plasma	/	CE 0123
W31-C	PSA	Cassette	4.0 mm	Serum/Plasma	/	CE 0123
W030-C	FOB	Cassette	4.0 mm	Feces	50 ng/mL	CE 0123 CLIA WAIVED
W30-S	FOB	Strip	2.5/3.0/4.0 mm	Feces	200 ng/mL	CE
W30-C	FOB	Cassette	4.0 mm	Feces	200 ng/mL	CE

Drug of Abuse Tests (Urine)

Cat. No.	Product	Format	Width	Specimen	Sensitivity	Qualification
W7-S	AMP (Amphetamine)	Strip	3.0/4.0 mm	urine	1000/500/300 ng/mL	CE CLIA WAIVED
W7-C		Cassette	4.0 mm	urine	1000/500/300 ng/mL	CE CLIA WAIVED
W7-P		Panel	4.0 mm	urine	1000/500/300 ng/mL	CE CLIA WAIVED
W8-S	BAR (Barbiturates)	Strip	2.5/3.0/3.5/4.0 mm	urine	300/200 ng/mL	CE CLIA WAIVED
W8-C		Cassette	4.0 mm	urine	300/200 ng/mL	CE CLIA WAIVED
W8-P		Panel	4.0 mm	urine	300/200 ng/mL	CE CLIA WAIVED
W9-S	BZO (Benzodiazepines)	Strip	2.5/3.0/3.5/4.0 mm	urine	300/200/100 ng/mL	CE CLIA WAIVED
W9-C		Cassette	4.0 mm	urine	300/200/100 ng/mL	CE CLIA WAIVED
W9-P		Panel	4.0 mm	urine	300/200/100 ng/mL	CE CLIA WAIVED
W10-S	COC (Cocaine)	Strip	2.5/3.0/3.5/4.0 mm	urine	300/150/100 ng/mL	CE CLIA WAIVED
W10-C		Cassette	4.0 mm	urine	300/150/100 ng/mL	CE CLIA WAIVED
W10-P		Panel	4.0 mm	urine	300/150/100 ng/mL	CE CLIA WAIVED
W11-S	mAMP/MET (Methamphetamine)	Strip	2.5/3.0/3.5/4.0 mm	urine	1000/500 ng/mL	CE CLIA WAIVED
W11-C		Cassette	4.0 mm	urine	1000/500 ng/mL	CE CLIA WAIVED
W11-P		Panel	4.0 mm	urine	1000/500/300 ng/mL	CE CLIA WAIVED
W12-S	MDMA (MDMA)	Strip	2.5/3.0/3.5/4.0 mm	urine	1000/500 ng/mL	CE CLIA WAIVED
W12-C		Cassette	4.0 mm	urine	1000/500 ng/mL	CE CLIA WAIVED
W12-P		Panel	4.0 mm	urine	1000/500 ng/mL	CE CLIA WAIVED
W13-S	BUP (Buprenorphine)	Strip	2.5/3.0/3.5/4.0 mm	urine	10 ng/mL	CE CLIA WAIVED
W13-C		Cassette	4.0 mm	urine	10 ng/mL	CE CLIA WAIVED
W13-P		Panel	4.0 mm	urine	10 ng/mL	CE CLIA WAIVED
W14-S	MOP (Morphine)	Strip	2.5/3.0/3.5/4.0 mm	urine	500/300/100 ng/mL	CE CLIA WAIVED
W14-C		Cassette	4.0 mm	urine	500/300/100 ng/mL	CE CLIA WAIVED
W14-P		Panel	4.0 mm	urine	500/300/100 ng/mL	CE CLIA WAIVED
W15-S	MTD (Methadone)	Strip	2.5/3.0/3.5/4.0 mm	urine	300/200 ng/mL	CE CLIA WAIVED
W15-C		Cassette	4.0 mm	urine	300/200 ng/mL	CE CLIA WAIVED
W15-P		Panel	4.0 mm	urine	300/200 ng/mL	CE CLIA WAIVED

Drug of Abuse Tests (Urine)

Cat. No.	Product	Format	Width	Specimen	Sensitivity	Qualification
W16-S	OPI (Opiate)	Strip	2.5/3.0/3.5/4.0 mm	urine	2000/300/100 ng/mL	CE CLIA WAIVED
W16-C		Cassette	4.0 mm	urine	2000/300/100 ng/mL	CE CLIA WAIVED
W16-P		Panel	4.0 mm	urine	2000/300/100 ng/mL	CE CLIA WAIVED
W17-S	PCP (Phencyclidine)	Strip	2.5/3.0/3.5/4.0 mm	urine	25 ng/mL	CE CLIA WAIVED
W17-C		Cassette	4.0 mm	urine	25 ng/mL	CE CLIA WAIVED
W17-P		Panel	4.0 mm	urine	25 ng/mL	CE CLIA WAIVED
W18-S	TCA (Tricyclic Antidepressants)	Strip	2.5/3.0/3.5/4.0 mm	urine	1000 ng/mL	CE CLIA WAIVED
W18-C		Cassette	4.0 mm	urine	1000 ng/mL	CE CLIA WAIVED
W18-P		Panel	4.0 mm	urine	1000 ng/mL	CE CLIA WAIVED
W19-S	THC (Marijuana)	Strip	2.5/3.0/3.5/4.0 mm	urine	50/40/25 ng/mL	CE CLIA WAIVED
W19-C		Cassette	4.0 mm	urine	50/40/25 ng/mL	CE CLIA WAIVED
W19-P		Panel	4.0 mm	urine	50/40/25 ng/mL	CE CLIA WAIVED
W21-S	OXY (Oxycodone)	Strip	2.5/3.0/3.5/4.0 mm	urine	100 ng/mL	CE CLIA WAIVED
W21-C		Cassette	4.0 mm	urine	100 ng/mL	CE CLIA WAIVED
W21-P		Panel	4.0 mm	urine	100 ng/mL	CE CLIA WAIVED
W55-S	KET (Ketamine)	Strip	2.5/3.0/3.5/4.0 mm	urine	1000/300 ng/mL	CE
W55-C		Cassette	4.0 mm	urine	1000/300 ng/mL	CE
W55-P		Panel	4.0 mm	urine	1000/300 ng/mL	CE
W57-S	PPX (Propoxyphene)	Strip	2.5/3.0/3.5/4.0 mm	urine	300 ng/mL	CE CLIA WAIVED
W57-C		Cassette	4.0 mm	urine	300 ng/mL	CE CLIA WAIVED
W57-P		Panel	4.0 mm	urine	300 ng/mL	CE CLIA WAIVED
W074-S	EDDP (EDDP)	Strip	2.5/3.0/3.5/4.0 mm	urine	300/100 ng/mL	CE
W074-C		Cassette	4.0 mm	urine	300/100 ng/mL	CE
W074-P		Panel	4.0 mm	urine	300/100 ng/mL	CE
W086-S	TRA (Tramadol)	Strip	2.5/3.0/3.5/4.0 mm	urine	1000/200 ng/mL	CE
W086-C		Cassette	4.0 mm	urine	1000/200 ng/mL	CE
W086-P		Panel	4.0 mm	urine	1000/200 ng/mL	CE
W087-S	FTY (Fentanyl)	Strip	2.5/3.0/3.5/4.0 mm	urine	20 ng/mL	CE
W087-C		Cassette	4.0 mm	urine	20 ng/mL	CE
W087-P		Panel	4.0 mm	urine	20 ng/mL	CE
W089-S	K2 (Synthetic Cannabis)	Strip	2.5/3.0/3.5/4.0 mm	urine	50 ng/mL	CE
W089-C		Cassette	4.0 mm	urine	50 ng/mL	CE
W089-P		Panel	4.0 mm	urine	50 ng/mL	CE
W139-S	COT (Cotinine)	Strip	2.5/3.0/3.5/4.0 mm	urine	100 ng/mL	CE
W139-C		Cassette	4.0 mm	urine	100 ng/mL	CE
W139-P		Panel	4.0 mm	urine	100 ng/mL	CE
W140-S	ETG (Ethyl Glucuronide)	Strip	2.5/3.0/3.5/4.0 mm	urine	500 ng/mL	CE
W140-C		Cassette	4.0 mm	urine	500 ng/mL	CE
W140-P		Panel	4.0 mm	urine	500 ng/mL	CE

Drug of Abuse Tests (Urine)

Cat. No.	Product	Format	Width	Specimen	Sensitivity	Qualification
W083-S	CAF (Caffeine)	Strip	2.5/3.0/3.5/4.0 mm	Urine	Variable	
W083-C		Cassette	4.0 mm	Urine	Variable	
W083-P		Panel	4.0 mm	Urine	Variable	
W155-S	MQL(Methaqualone)	Strip	4.0 mm	Urine	Variable	
W155-P		Panel	4.0 mm	Urine	Variable	

Multi-Drug Tests

Cat. No.	Product	Format	Width	Test items	Specimen	Sensitivity	Qualification
W6506	Oral Fluid Multi-Drug 6 Test	Cassette	/	6 Items (AMP/BZO/COC/THC/MET/MOP)	Oral Fluid	Variable	/
W602~610-CU6	Multi-Drug T-Cube Test	T-Cube	3.0 mm	2-10 Items	Oral Fluid	Variable	/
W200X(2-12)-P	Urine Multi-Drug Test	Panel	4.0 mm	2-12 Items	Urine	Variable	☹☹ CLIA WAIVED
W200X(13-16)-P		Panel	4.0 mm	13-16 Items	Urine	Variable	☹☹
W50X(2-18)-CU2		T-Cup	2.5 mm	2-18 Items	Urine	Variable	☹☹
W50X(2-18)-CU3		Q-Cup	2.5 mm	2-18 Items	Urine	Variable	☹☹
W50X(2-18)-CU8		T-Cup Plus	2.5 mm	2-18 Items	Urine	Variable	☹☹
W50X(2-18)-CU7		Q-Cup Plus	2.5 mm	2-18 Items	Urine	Variable	☹☹

Other Rapid Tests

Cat. No.	Product	Format	Width	Specimen	Sensitivity	Qualification
W5410-S	Urinalysis Reagents Strip (1~10 Parameter)	Strip	5.0 mm	Urine	/	☹☹
W46-C4P	Troponin I Test	Cassette	4.0 mm	S/P/WB	1 ng/mL	☹☹
W53-S	Alcohol Saliva Test	Strip	3.5/4.0 mm	Oral Fluid	Custom	☹☹

Dry Chemistry Analyzer DC-101

- Liver Function
- Kidney Function
- Emergency Diagnosis
- Myocardial Enzyme
- Blood Routine Examination
- Blood Glucose



ACCURATE

Advanced reflective photometry method ensures accurate result, good correlation with reference method.

FLEXIBLE

- Small & Portable
- Required less sample volume ≤ 20 µL

RAPID

- 60 tests per hour
- Result comes out in 1-5 min(s)

EASY-TO-USE

- One step loading
- Full-touch screen, user-friendly interactive system.

ECONOMY

- Single test, free combination, no waste.
- No pipeline system, maintenance free.

Test Item

Catalog No.	Category	Test Item	Sample Type	Reaction Time(min)	Sample Volume(µL)	Measuring Range
W416	Liver Function	ALT	Serum/Plasma	2.5	15	5-800 U/L
W420		AST	Serum	2	15	10-800 U/L
W414		ALB	Serum	1	15	10-60 g/L
W412	Kidney Function	UREA	Serum/Plasma	2	15	2.5-40 mmol/L
W418		CREA	Serum	5	15	30-1500 µmol/L
W413		UA	Serum/Plasma	2	15	120-1200 µmol/L
W422	Myocardial Enzyme	CK	Serum	3	15	20-1600 U/L
W423		CK-MB	Serum	4	15	5-300 U/L
W419		LDH	Serum	2	15	100-1800 U/L
W426		a-HBDH	Serum	2	15	20-1000 U/L
W411	Emergency Treatment	α-AMY	Serum/Urine	2	10	20-2000 U/L
W424		LPS	Serum	3	10	5-900 U/L
W417		CHE	Serum	2	15	200-12500 U/L
W421		ALC	Serum	2	15	10.0-300.0 mg/dL
W415	Blood Routine Examination	HB	Whole Blood	1	10	45-200 g/L
W425	Blood Glucose	GLU	Serum	2	15	1.00-25.00 mmol/L

* Ongoing: TP, TBIL, DBIL, TBA, G6PD;
TC, TG, LDL-C, HDL-C, DAO, D-LAC, BET, NH3

Blood Gas Analyzer BGA-101



pH, pCO₂, pO₂, K⁺, Na⁺, Ca⁺⁺, Cl⁻, Hct, Glu, Lac

Calculation Parameters: cH⁺, cH⁺(T), pH(T), pCO₂(T), pO₂(T), HCO₃⁻atd, HCO₃⁻std, BB(B), BE(B), BE(ecf), ct CO₂, Ca⁺⁺(7.4), AnGap, tHb(est), sO₂(est), pO₂(A-a), pO₂(A-a)(T), pO₂(a/A), pO₂(a/A)(T), RI, RI(T), pO₂/FI O₂, pO₂(T)/FI O₂

Specifications:

- Sample Type: Whole Blood
- Sample Volume: 80 µL
- Reading Time: 30 s
- Suitable for Syringe and Capillary
- Memory:>50,000
- Good Correlation to Reference Method: R²>0.95
- Test card is stable and effective for 15 months
- Reagent pack can be used for up to 60 days after installation



Test Packages

Test Item \ Cat No.	W459	W460	W461	W462	W463	W464	W465	W466
pH	•			•	•	•	•	•
pCO ₂	•			•	•	•	•	•
pO ₂	•			•	•	•	•	•
K ⁺		•	•	•	•	•	•	•
Na ⁺		•	•	•	•	•	•	•
Ca ⁺⁺			•			•	•	•
Cl ⁻		•	•	•	•		•	•
Hct					•	•	•	•
Glu								•
Lac								•

Optical Coagulation Analyzer OCG-102

Cat. No.	Item	Specimen	Specification	Qualification
OCG-102	Optical Coagulation Analyzer	/	/	CE
W480-484	PT(INR)/APTT/TT/FIB/ACT	Citrated Venous Whole Blood and Capillary Blood (PT/INR only)	24 T/Kit	CE

CE

PT/INR, APTT, TT, FIB, ACT

Specifications:

- **Sample Type:** Citrated Venous Whole Blood or Capillary Blood (PT/INR Only)
- **Sample Volume:** 20 µL
- **Time to Result:** ≤5 minutes
- **Precision:** CV≤5%
- **Store at room temperature and valid for 18 months**



Blood Glucose Monitoring System

Cat. No.	Item	Specimen	Sample Volume
WM801W02	Blood Glucose Monitoring System	/	/
W301	Blood Glucose Test Strip	Capillary Whole Blood	8 µL

• **Glucose Meter**



Immunofluorescence Quantitative Test System

• Finecare™ Platform



Benefits and Advantages:

- Reliable
- Fast and Accurate, Result in 3-15 mins, Full Traceability, QC Available
- Cost-efficient
- Full Test Line, CE/CFDA/FSC/NGSP/IFCC Certified, Competitive Price
- Easy-to-use
- LCD Touch Screen, Multi-lingual, Internal Thermal Printer, Built-in Temperature Control
- Room Temperature Stable Test Kit
- Store at 2-30 C, 24 Months Expiry Period
- Date Management and Transmission
- HIS/LIS Connectivity

Finecare™ FIA Meter serise is a fluorescence immunochromatographic analysing system which can help diagnose conditions such as infection, diabetes, cardiovascular diseases, renal injury and cancers, etc.

The Finecare™ platform provides reliable and quantitative results of various kinds of analytes in human blood or urine within several minutes.


Application: Laboratory, ER, Cardiology, ICU, Respiratory, Pediatrics, etc.

Cat. No.	Category	Test Item	Diseases	Specimen Type	Sample Volume	Reaction Time	Measuring Range	Qualification
W276	COVID-19	SARS-CoV-2 Ab ^{New}	COVID-19	S/P/WB	10 µL	10 min	Qualitative	CE
W277		SARS-CoV-2 IgM ^{New}	COVID-19	S/P/WB	10 µL	10 min	Qualitative	CE
W241	Vitamin	Vitamin D	Bone Metabolism	Plasma/Serum	75 µL	15 min	5-100 ng/mL	CE
W207	Diabetic & Renal Injury Markers	HbA1c	Diabetes Mellitus	Whole Blood	10 µL	5 min	4.0-14.5%	CE IFCC NGSP
W206		MAU	Renal Injury	Urine	75 µL	3 min	5-300 mg/L	CE
W219		CysC	Renal Injury	S/P/WB	10 µL	5 min	0.2-10 mg/L	CE
W228		NGAL	Renal Injury	Urine	75 µL	15 min	10-1500 ng/mL	CE
W229		β2-MG	Proximal Renal Tubule Dysfunction	S/P/WB	10 µL	10 min	0.3-20 mg/L	CE
W231	Thyroid Function	T3	Thyroid Malfunction	S/P/WB	75 µL	15 min	0.61-9.22 nmol/L	CE
W232		T4	Thyroid Malfunction	S/P/WB	75 µL	15 min	12.87-300 nmol/L	CE
W220		TSH	Thyroid Malfunction	S/P/WB	75 µL	15 min	0.1-100 mIU/L	CE
W200		ft3 ^{New}	Thyroid Malfunction	S/P/WB	75 µL	15 min	0.40-50.00 pmol/L or 0.26~32.55 ng/dL	CE
W259		ft4 ^{New}	Thyroid Malfunction	S/P/WB	75 µL	15 min	1.00~100 pmol/L or 0.078~7.77 ng/dL	CE
W225	Fertility	β-HCG	Fertility	S/P/WB	20 µL	15 min	2-200,000 mIU/mL	CE
W244		LH	Fertility	S/P/WB	75 µL	15 min	1-100 mIU/mL	CE
W245		FSH	Fertility	S/P/WB	75 µL	15 min	1-100 mIU/mL	CE
W233		Prog	Fertility	S/P/WB	75 µL	15 min	1.5-60 ng/mL	CE
W246		PRL	Fertility	S/P/WB	75 µL	15 min	1-200 ng/mL	CE
W247		E2	Fertility	S/P/WB	75 µL	15 min	9-3000 pg/mL	CE
W248		Testosterone	Fertility	S/P/WB	75 µL	15 min	0.2-15 ng/mL	CE
W866		AMH ^{New}	Fertility	S/P/WB	75 µL	15 min	0.1-16 ng/mL	CE
W202	Cardiac Markers	NT-proBNP	Heart Failure	S/P/WB	75 µL	15 min	18-35000 pg/mL	CE
W203		cTn I	Myocardial Infarction	S/P/WB	75 µL	15 min	0.1-50 ng/mL	CE
W204		Myo	Myocardial Infarction	S/P/WB	75 µL	15 min	2.0-400 ng/mL	CE
W205		CK-MB	Myocardial Infarction	S/P/WB	75 µL	15 min	0.3-100 ng/mL	CE
W217		H-FABP	Myocardial Infarction	S/P/WB	75 µL	15 min	1-120 ng/mL	CE
W216		3 in 1 (cTn I/Myo/CK-MB)	Myocardial Infarction	S/P/WB	75 µL	15 min	Same with single items	CE
W222		BNP	Heart Failure	S/P/WB	75 µL	15 min	5-5000 pg/mL	CE
W223		2 in 1 (cTn I/ NT-proBNP)	Myocardial Infarction	S/P/WB	75 µL	15 min	Same with single items	CE
W240		hsCRP	CVD	S/P/WB	10 µL	3 min	0.1mg/L-10mg/L	CE
W254		TnT ^{New}	Myocardial Infarction	S/P/WB	100 µL for whole blood 75 µL for plasma/serum	15 min	0.03-10 ng/mL	CE
W242	Lp-PLA2	Atherosclerosis	S/P/WB	75 µL	15 min	10-1000 ng/mL	CE	
W239	Craniocerebral Injury Markers	S100β	Cerebral Ischemia Cerebral Hemorrhage	S/P/WB	75 µL	15 min	0.05-10 ng/mL	CE
W211	Coagulation	D-Dimer	Thrombosis	WB/Plasma	75 µL	5 min	0.1-10 mg/L	CE
W227		One step D-Dimer	Thrombosis	WB/Plasma	100 µL for whole blood 75 µL for plasma/serum	15 min	0.1-10 mg/L	CE
W210	Inflammation Markers	PCT	Sepsis	S/P/WB	50/75 µL	15 min	0.1-100 ng/mL	CE
W201		CRP (hsCRP+CRP)	Inflammation	S/P/WB	5/8.5 µL	3 min	0.5-200 mg/L	CE
W218		2 in 1 (CRP+PCT)	Infection & Inflammation	S/P/WB	75 µL	15 min	CRP:0.5-150 mg/L PCT:0.1-100 ng/L	CE
W251		SAA	Infection & Inflammation	S/P/WB	10 µL	15 min	1.0-300 mg/L	CE
W221		IL-6	Inflammation	S/P/WB	75 µL	15 min	3-4000 pg/mL	CE
W208	Tumor Markers	AFP	Liver Cancer	S/P/WB	75 µL	15 min	5-400 ng/mL	CE
W209		PSA	Prostate Cancer	S/P/WB	75 µL	15 min	2-100 ng/mL	CE 0123
W226		CEA	Cancer	S/P/WB	75 µL	15 min	1-500 ng/mL	CE
W230		fPSA	Prostate Cancer	S/P/WB	75 µL	15 min	0.2-30 ng/mL	CE 0123
W250	Others	Cortisol	Cortisol disorders	S/P/WB	75 µL	15 min	50 nmol/L~1000 nmol/L	CE

Wondfo®

OUR VISION

Provide professional products and services of rapid diagnosis and chronic disease management; Pursue better and convenient user experience; Promote high-standard life and health quality; Fulfill shared dream of our staff, customers and partners.



1,500+ Employees

200+ R&D Staffs

2 Overseas Lab

1 Branch

10+ Overseas Offices

WONDFO BIOTECH
WeAreWorkingForYourHealth

Wondfo®

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Luogang District,
510663, Guangzhou, P.R. China

Tel: +86-20-32296083
Fax: +86-20-32296063
E-mail: sales@wondfo.com.cn
Website: www.en.wondfo.com.cn

Customers Service Center:
Tel: +86-20-32296083
Fax: +86-20-32296063
E-mail: customerservice@wondfo.com.cn



物料编码及项目名称：41000944 国际产品总目录（2020）



尺寸规格：210*285mm

颜色：四色

设计师：许广添

材质：封面157克双铜过哑膜内页80克双铜+骑马钉 申请人：

设计时间：2020.03.27
