



# Finecare™ FIA METER III PLUS

## Multi-Channel

20 channels, detect more items one time

## Reliable Accuracy

Good correlation with CLIA method

## Temperature control

Provide better reaction condition

## Rapid Test

High speed of detection, enable 200 tests per hour

## FS-205

- Dimensions: 374\*396\*480.5 mm
- Weight: Approx. 19.5Kg

## ID chip inlet

- Automatically reading of batch number and calibration information in the ID chip

## Test Channel

- 20 channels
- Random access for samples
- Automatic internal incubation

## Built-in Thermal Printer

- Real-time auto-printing
- Support for external printers

## Designed Sampler

- Specially-designed sampler for small volume of specimen

## Large Display

- 10" LCD Touch Screen
- Intuitive user interface

## Connection

- Direct connection to LIS/HIS


## Cartridge Tray

- Cartridges will be automatically ejected into the tray when test is finished


## Test Reagent

- Room temperature storage
- Suitable for whole blood, serum, plasma and urine
- Get results in 3-15 minutes


## TEST PROCEDURE




Sampling



Mixing



Dropping



Reading

## APPLICATION



Laboratory



ER



ICU



Clinics



GP Office



物料编码及项目名称：41000896 飞测IIIPlus双折页2020

申请人：李朦

尺寸规格：420\*285mm

颜色：四色

设计师：杨晓洁

材质：200克哑粉纸

设计时间：2020.03.12