

HIV Self-Test



- **Less blood**

Just needs 10 μ L of blood from your fingertip.

- **Privacy protection**

Get results from the privacy of your own home.


- **Easy-to-use**

Test in 3 steps
Receive result in 15 mins.




OPERATION PROCEDURE


Step 1




Step 2




Step 3



Result Interpretation



15~30 min



Clinical Performance

Users	Sensitivity	Specificity
Professional users	100%	99.96%
Self-test users	95.8%	99.6%

Order Information

Catalog No.	Pack size
W006P0058	1 T/kit
W006P0059	20 T/kit
W006P0060	100 T/kit

Guangzhou Wondfo Biotech Co., Ltd.

Add: No. 8 Lizhishan Road, Science City, Luogang District, 510663, Guangzhou, P.R. China
 Tel: +86-20-32296082 Fax: +86-20-32296062
 E-mail: global@wondfo.com.cn Websites: en.wondfo.com

