

# SAMPLE COLLECTION KIT

Various swab configurations of 1, 5, 10 and 20 applicators with one tube 3 mL VTM are available for individual or pooled sample requirement from different clinical settings.



## Safety

Good sealing performance of sampling tube ensures no liquid leakage.



## Efficient

Optimized virus preservation solution and flocked swab ensures the effective release of virus from clinical specimen.



## Compatible

A compatible collection and preservation kit of various viruses, such as 2019-nCoV and influenza.



## VTM Tube

The vial is made of polypropylene allows for safe collection and transport of specimen.



## VTM Solution

VTM solution contains guanidine composition to inactivate virus effectively.



## Swab

High-quality flocked swabs. 3 cm and 8 cm moded breakpoints compatible with oropharyngeal and nasopharyngeal swab.

## ORDER INFORMATION

Cat. No	Product Type	Package Size	Storage Condition	Shelf Life
WH01P0007	Inactivated VTM with Nasopharyngeal swab	50 T/Box	2~30°C	24 M
WH01P0008	Non-Inactivated VTM with Nasopharyngeal swab			
WH01P0009	Inactivated VTM with Oropharyngeal swab			
WH01P0010	Non-Inactivated VTM with Oropharyngeal swab			

