

One Step FSH Urine Test

Easy way to determine the menopause status!

Wondfo One Step FSH Urine Test is a qualitative test used to detect the levels of FSH associated to the symptoms and stages of menopause. It can help you to determine if your hot flashes, irregular periods and other health symptoms are menopause related. High FSH level in women may indicate menopause or infertility.



RELIABLE

High-quality service is provided.

CONVENIENT

Easy to use.

RAPID


Fast result within 10 minutes.

WARM PROMPT:


LH and FSH work together to promote the follicular development and ovulation. You can use Wondfo One Step Ovulation Urine Test to double check as well.

TEST PROCEDURE & INTERPRETATION


Strip



Cassette



Midstream



ORDER INFORMATION

Catalog No.	Format	Width	Specimen	Sensitivity	Test time	Storage temperature	Shelf life	Qualification
W52-S	Strip	3.5mm	Urine	25 mIU/mL	10 min	4-30°C	24 months	CE 0123
W52-C	Cassette	4.0mm						
W52-M	Midstream	6.0mm						

