

HCG One Step HCG Urine Test

Wondfo One Step HCG Urine Test is a rapid, qualitative detection of human chorionic gonadotropin ("hCG") which appears in urine very early during pregnancy.



RELIABLE

High-quality service is provided.

CONVENIENT


Easy to use.

RAPID

Fast result within 3 minutes.

TEST PROCEDURE & INTERPRETATION


Strip



Cassette



Midstream



ORDER INFORMATION

| Catalog No. | Format | Width | Specimen | Sensitivity | Test time | Storage temperature | Shelf life | Qualification |
|---|-----------|---------------|----------|--------------------|-----------|---------------------|------------|---|
| W1-S/ W1-S(20mIU)/ W1-S(10mIU) | Strip | 3.5mm | Urine | 25/20/10 mIU/mL | 3 min | 4-30°C | 36 months | FDA (Excluding 20 mIU/mL) |
| W1-C/ W1-C/ W1-C(20mIU)/ W1-C(10mIU) W1-C(10mIU)/ W1-C(10mIU) | Cassette | 3.5/4.0mm | | | | | | |
| W1-M/ W1-MII/ W1-MIII/ W1-M4/ W1-M6/ W1-M(20mIU)/ W1-MII(20mIU)/ W1-MIII(20mIU)/ W1-M4(20mIU)/ W1-M6(20mIU)/ W1-M(10mIU)/ W1-MII(10mIU)/ W1-MIII(10mIU)/ W1-M4(10mIU)/ W1-M6(10mIU) | Midstream | 3.5/5.0/6.0mm | | | | | | |



ONE STEP HCG URINE/SERUM TEST

Professional way to determine the pregnancy status!

The Wondfo One Step HCG Urine/Serum Test in vitro diagnostic test is used for the qualitative determination of hCG in urine or serum for early detection of pregnancy. It detects the presence of human chorionic gonadotropin (“hCG”) at a concentration of 25 mIU/mL or above.

For professional use only.



RELIABLE

High-quality service is provided.

CONVENIENT

Easy to use.

RAPID

Fast result within 5 minutes.


ACCURATE

Over 99% accuracy.


TEST PROCEDURE

Strip

* Do not allow the urine level to exceed the MAX line.



Cassette



ORDER INFORMATION

| Catalog No. | Packing | Format | Width | Specimen | Sensitivity | Test time | Storage temperature | Shelf life | Qualification |
|----------------------|--------------|----------|-------|--------------|-------------|-----------------------------|---------------------|------------|---------------|
| W1-S3P/ W001P0720 | 50/ 100T/Box | Strip | 3.5mm | Urine/ Serum | 25 mIU/mL | 3min(Urine)/ 5min(Serum) | 4-30°C | 24 months | CE FDA |
| W1-C3P/ W001P0721 | 25/ 40T/Box | Cassette | 4.0mm | Urine/ Serum | 25 mIU/mL | 3min(Urine)/ 5min(Serum) | 4-30°C | 24 months | CE FDA |

Guangzhou Wondfo Biotech Co., Ltd.

No.8 Lizhishan Road, Science City, Luogang District, 510663, Guangzhou, P.R.China

Tel: (+86)400-830-8768

E-mail: sales@wondfo.com.cn

Website: en.wondfo.com.cn

Wondfo reserves all rights. All contents of this brochure are copyrighted by Wondfo.

