



WONDFO BIOTECH
WeAreWorkingForYour **Health**

HbA1c

(Hemoglobin A1c)

Finecare[™]

HbA1c RAPID QUANTITATIVE TEST
(Diabetes Markers)

For Better and Easier Glycemic Control



Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City,
Luogang District, 510663 Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com.cn
E-mail: sales@wondfo.com.cn




Rev. A2 Rel.: 2021/12/15

Finecare™

HbA1c Rapid Quantitative Test

- ✓ Help diagnosis and monitoring of diabetes
- ✓ Just 10 µL blood required from fingerstick or venous blood
- ✓ Room temperature storage up to 24 months
- ✓ Certified by **CE**, **NGSP**, **IFCC**, **RIQAS**
- ✓ Traceable to the IFCC reference system
- ✓ Selectable results reporting (% DCCT, IFCC mmol/mol, eAG mg/dl) in 5 minutes




CLINICAL APPLICATION

- Diabetes or prediabetes screening and diagnosis
- Helps monitoring and chronic management of diabetes
- Helps gestational diabetes control
- Helps prescribing and monitoring the medicines for diabetes
- Improves clinical decision-making at primary and point of care

ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Loading Volume	Reaction Time	Detection Range	Storage
HbA1c	W207	Fingerstick Blood / Venous Whole Blood	10 µL	75 µL	5 minutes	4~14.5%	4~30°C for 24 months

TEST PROCEDURES



ADAPTATION INSTRUMENTS

1. Finecare™ FIA Meter (FS-112)
2. Finecare™ FIA Meter Plus (FS-113)
3. Finecare™ II Plus SE (FS-114)
4. Finecare™ III Plus (FS-205)



Point of Care
(GP, Doctor's Office, Clinics,
Pharmacy, Family Physician)



Laboratory



In-patient Care



Hospital Clinic
Department



Diabetes Clinic or
Diabetes Center