



WONDFO BIOTECH
WeAreWorkingForYour **Health**

iFOB

FinecareTM

iFOB RAPID QUANTITATIVE TEST

Instrument detection, automatic display of quantitative results



Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City,
Luogang District, 510663 Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com.cn
E-mail: sales@wondfo.com.cn




Rev. A1 Rel.: 2022/01/11

Finicare™

iFOB Rapid Quantitative Test

- ✓ Help detecting gastrointestinal bleeding
- ✓ Help in the initial screening of colorectal cancer
- ✓ Convenient all-in-one sampler included
- ✓ Results are ready in only 5 minutes
- ✓ Room temperature storage up to 24 months



CLINICAL APPLICATION

Fecal Occult Blood Test (FOBT) is mainly used for initial screening of colorectal cancer.


The ACS (American Cancer Society) recommends that people at average risk of colorectal cancer start regular screening at age 45.

Cancer Council recommends all Australians aged 50 and over screen for bowel cancer with an FOBT every two years.

ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sampling	Loading Volume	Reaction Time	Detection Range	Storage
iFOB	W279P0001	Feces	6 different parts of the feces	3 drops	5 minutes	25~1000 ng/mL	4~30°C for 24 months

TEST PROCEDURE



ADAPTATION INSTRUMENTS

1. Finicare™ FIA Meter Plus (FS-113)
2. Finicare™ II Plus SE (FS-114)
3. Finicare™ III Plus (FS-205)



Point of Care
(GP, Doctor's Office, Clinics,
Pharmacy, Family Physician)



Laboratory



Anorectal
Department



Gastroenterology
Department



Physical
examination
Department

