



WONDFO BIOTECH
WeAreWorkingForYour **Health**

PCT

(Procalcitonin)

Finecare[™]

PCT RAPID QUANTITATIVE TEST

(Infection Markers)

For Timely Sepsis Diagnosis and Better Antibiotics Use

Wondfo

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City,
Luogang District, 510663 Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com.cn
E-mail: sales@wondfo.com.cn




Rev. A1 Rel.: 2021/12/17

Finicare™


PCT Rapid Quantitative Test

- ✓ High sensitivity and specificity for bacterial infection
- ✓ The biomarker of choice to early and effective diagnosis of sepsis
- ✓ Guiding the use of antibiotics in primary care and critical care
- ✓ Room temperature storage up to 24 months
- ✓ Results are ready in only 15 minutes



CLINICAL APPLICATION

- Early diagnosis of severe infection
- Differentiating infection of bacteria from virus



- Accurately guiding medication time and dosage of antibiotics
- Effectively avoiding antibiotic abuse


- Severity monitoring of severe patient
- Monitoring of postoperative infection

- Sepsis indication, MODS prognosis
- Clinical observation of antibiotic treatment

ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Loading Volume	Reaction Time	Detertion Range	Storage
PCT	W210 / W210P0004	Whole blood/ Serum/Plasma	75 µL of whole blood or 50 µL serum or plasma	75 µL	15 minutes	0.1-100 ng/mL	4-30°C for 24 months

TEST PROCEDURE



ADAPTATION INSTRUMENTS

1. Finicare™ FIA Meter Plus (FS-113)
2. Finicare™ II Plus SE (FS-114)
3. Finicare™ III Plus (FS-205)

