

Wondfo



WONDFO BIOTECH
WeAreWorkingForYour **Health**

Thyroid Function

Finecare[™]

Thyroid Function Rapid Quantitative Test

Regular screening of thyroid function!

Wondfo

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City,
Luogang District, 510663 Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com.cn
E-mail: sales@wondfo.com.cn



Rev. A1 Rel.: 2022/08/15

Finecare™

Thyroid Function Rapid Quantitative Test

- ✓ Direct detection of Thyroid hormone levels
- ✓ Help diagnosis Graves' Disease, Hashimoto's Disease, Congenital hypothyroidism (CH)
- ✓ Support whole blood samples; Results ready in 15 minutes
- ✓ Room temperature storage up to 24 months

CLINICAL APPLICATION

- Diagnose and monitor thyroid diseases
- Find the cause of an underactive thyroid gland
- Keep track of thyroid function in hyperthyroidism treatment
- Double-check the diagnosis of congenital hypothyroidism

APPLICATION SCENARIOS

Point of Care (GP, Doctor's Office, Clinics) Lab Endocrinology Department Gynaecology and Obstetrics Chest pain center

ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Reaction Time	Detection Range	Storage
T3	W231	Whole blood/ Serum/Plasma	75 µL	15 minutes	0.61~9.22 nmol/L (0.4~6.0 ng/mL)	4~30°C for 24 months
T4	W232				12.87~300 nmol/L	
TSH	W220				0.1~100 mIU/L	
fT3	W200				0.40~50.00 pmol/L or 0.26~32.55 pg/mL	
fT4	W259				1.00~100 pmol/L or 0.078~7.77 ng/dL	

TEST PROCEDURE

ADAPTATION INSTRUMENTS

1. Finecare™ FIA Meter Plus (FS-113)
2. Finecare™ II Plus SE (FS-114)
3. Finecare™ III Plus (FS-205)

