

MALARIA

One Step Malaria Test

Malaria is a lethal parasite disease, spread through mosquito. There are 4 species of malaria species invading human being: PV, PF, PM and PO. They collectively referred to as Pan.



Guangzhou Wondfo Biotech Co., Ltd.

No.8 Lizhishan Road, Science City, Luogang District,

510663, Guangzhou, P.R.China

E-mail: sales@wondfo.com.cn

Tel: (+86) 400-830-8768

Website: en.wondfo.com




WONDFO MALARIA TEST

Product Name	Format	Pack Size	Qualification
One Step Malaria HRP2 (P.f) Test	Cassette	25 Tests/kit; 40 Tests/kit	CE
	Strip	50 Tests/kit; 100 Tests/kit	CE
One Step Malaria HRP2/pLDH (P.f/P.v) Test	Cassette	25 Tests/kit; 40 Tests/kit	CE
One Step Malaria HRP2/pLDH (P.f/Pan) Test	Cassette	25 Tests/kit; 40 Tests/kit	CE

CHARACTERISTICS

Specimen <hr/> Venous whole blood/ Whole blood collected by fingerstick	Reaction Time <hr/> 15 mins	Storage Temperature <hr/> 2~35 °C	Shelf Life <hr/> 24 months
-----------------------------------------------------------------------------------------	----------------------------------------------	----------------------------------------------------	---------------------------------------------

TEST PROCEDURE



P.f	P.f / P.v		P.f / Pan	
Positive	P.f positive(+)	P.v positive(+)	P.f positive(+)	Non-P.f positive(+)
Negative	Negative(-)	Mixed infection with P.f/Pv(+)	Negative(-)	P.f monoinfection(+) or Mixed infection
Invalid	Invalid	Invalid	Invalid	Invalid
Invalid	Invalid	Invalid	Invalid	Invalid