


# One Step Syphilis Test

Wondfo One Step Syphilis Whole Blood/Serum/Plasma Test is a rapid test for the visual detection of antibodies of Treponema Pallidum (TP) as an aid in the diagnosis of syphilis infection.



## SENSITIVITY

100%

## SPECIFICITY

98%

## FORMAT

Cassette / Strip

## SHELF-LIFE

24 months

## SPECIMEN

Whole Blood/Serum/Plasma


## STORAGE TEMPERATURE

4°C~30°C




## TEST PROCEDURE

### For Cassette



### For Strip



## PERFORMANCE CHARACTERISTICS

Wondfo kit	TPPA		Total
	Positive	Negative	
Positive	53	1	54
Negative	0	50	50
Total	53	51	104

**Sensitivity of Wondfo One Step Syphilis Whole Blood/Serum/Plasma Test: 100%**

**Specificity of Wondfo One Step Syphilis Whole Blood/Serum/Plasma Test: 98%**

## ORDER INFORMATION

Item Code	Format	Description	Qualification
W34-C4P/W034P0071	Cassette	One Step Syphilis Whole Blood/Serum/Plasma Test	CE
W34-S4P/W034P0070	Strip	One Step Syphilis Whole Blood/Serum/Plasma Test	CE