

Category	Cat. No	Test Item	Specimen Type	Specimen Volume	Reaction Time	Quality Control	Qualification
Cardiac Markers	WS02P0001	One Step hs-cTnI/CK-MB/Myo	WB/S/P	75µL Serum/Plasma 100µL Whole Blood	8 mins	✓	CE
	WS01P0001	One Step hs-cTnI	WB/S/P	75µL Serum/Plasma 100µL Whole Blood	8 mins	✓	CE
	WS04P0001	One Step NT-proBNP	WB/S/P	75µL Serum/Plasma 100µL Whole Blood	8 mins	✓	CE
	WS03P0001	One Step BNP	WB/S/P	75µL Serum/Plasma 100µL Whole Blood	8 mins	✓	CE
	WS06P0001	One Step D-Dimer	WB/P	75µL Plasma 100µL Whole Blood	8 mins	✓	CE
Inflammation	WS08P0001 WS08P0002	CRP	WB/S/P	5µL Serum/Plasma 8.5µL Whole Blood	3 mins	✓	CE
	WS12P0001	One Step PCT	WB/S/P	75µL Serum/Plasma 100µL Whole Blood	3 mins	✓	CE
	WS14P0001	One Step IL-6	WB/S/P	75µL Serum/Plasma 100µL Whole Blood	3 mins	✓	CE
	WS09P0001 WS09P0002	SAA	WB/S/P	5µL Serum/Plasma 8.5µL Whole Blood	3 mins	✓	CE
Infectious Disease	Adeno		Ongoing				
	Adeno/RSV						
	Covid-19 Ag						
	Covid-19 Ag/FluA+B						
	FluA/B						
	Mycoplasma						
	RSV						
Strep A							

Category	Cat. No	Item	Format	Level	Qualification	Specification
Quality Control	WU02P0001	hs-cTnI/CK-MB/Myo Control	Lyophilized	3	CE	0.5mL*3 vials/box
	WU01P0001	hs-cTnI Control	Lyophilized	3	CE	0.5mL*3 vials/box
	WU04P0001	NT-proBNP Control	Lyophilized	3	CE	0.5mL*3 vials/box
	WU03P0001	BNP Control	Lyophilized	3	CE	0.5mL*3 vials/box
	WU05P0001	D-Dimer Control	Lyophilized	3	CE	0.5mL*3 vials/box
	WU06P0001	CRP Control	Lyophilized	3	CE	0.5mL*3 vials/box
	WU06P0002		Lyophilized	3	CE	0.2mL*3 vials/box
	WU09P0001	PCT Control	Lyophilized	3	CE	0.5mL*3 vials/box
	WU09P0002		Lyophilized	3	CE	0.2mL*3 vials/box
	WU10P0001	IL-6 Control	Lyophilized	3	CE	0.5mL*3 vials/box
	WU10P0002		Lyophilized	3	CE	0.2mL*3 vials/box
	WU07P0001	SAA Control	Lyophilized	3	CE	0.5mL*3 vials/box
	WU07P0002		Lyophilized	3	CE	0.2mL*3 vials/box



Finecare™




FS-2000

Fluorescence Immunoassay Analyzer



FS-2000 Fluorescence Immunoassay Analyzer >




Specification

- Dimension: 325*475*395 mm
- Weight: 19 kg



Connectivity

- WIFI connection to Bi-directional LIS/HIS
- RS232 and ethernet port
- 5 USB ports



Built-in Code Scanner

- No need for ID chip
- Sample batch information



Built-in Waste Bin

- Automatic storage of used cartridges
- Holds up to 100 used cartridges



Special cooling function

- Ensures stability under high temperature



Rechargeable Internal Battery

- 10400 mAh battery
- 24-hour standby
- Low battery alert



Large Display

- 10-inch capacity touch screen
- Intuitive user interface
- Smooth operation experience like smart-phone



Built-in Thermal Printer

- Real-time auto printing
- External printer supported



Linux System

- Smooth operation
- Security system
- Custom intelligent management module
- Multi-user management

Technical Advantages >



Stability

Cutting edge TRFIA methodology
High stability, lower CV



Convenient

Simple operation procedure
Automatic collection of used cartridges



Instant Results

Results delivered in less than
8 minutes



High Efficiency

Maximum simultaneous detection
of up to 16 items



Delicate Design

Automatic reagent door
Charging indicator light
Forward USB port

1 QR Code Technology


- No need manual input cartridge information
- Reagent card recognition success rate is as high as 99%

2 Smaller Cartridge

- Lower shipping costs

3 Excellent CV and Accuracy with Optimized FIA Methodology

- Satisfying and stable results delivered in short time



Operation Procedure >



Entry Test Information



Sample Loading



Start Testing

Application Scope >



Ambulance



Emergency Room



Clinics



Pharmacy



GP Office