

Finecare™ FIA METER PLUS

Convenient

Portable
Easy to operate

Reliable Accuracy

Good correlation
with CLIA method

Temperature control

Provide better
reaction condition

Rapid Test

Get results in
3-15 minutes

FS-113

- Dimension: 270*238*146mm
- Weight: 2.8 Kg

Connection

- RS 232 computer interface port/ Ethernet network

ID chip inlet

- Automatically reading of batch number and calibration information in the ID chip

Designed Sampler

- Specially-designed sampler for small volume of specimen

Inner Temperature Control

Built-in Thermal Printer

- Real-time auto-printing
- Support for external printers

Large Display

- 8" LCD Touch Screen
- Intuitive user interface


Test Channel

- Multiple testing modes available
- Unique quick mode for large number of samples


Test Reagent

- Room temperature storage
- Suitable for whole blood, serum, plasma and urine


TEST PROCEDURE




Sampling



Mixing



Dropping



Reading



Laboratory



ER



ICU



Clinics



GP Office

APPLICATION