

Finecare™ FIA METER III PLUS

Multi-Channel	Reliable Accuracy	Temperature control	Rapid Test
20 channels, detect more items one time	Good correlation with CLIA method	Provide better reaction condition	High speed of detection, enable 200 tests per hour

FS-205

- Dimensions:374*396*480.5 mm
- Weight: Approx. 19.5Kg

ID chip inlet

- Automatically reading of batch number and calibration information in the ID chip

Test Channel

- 20 channels
- Random access for samples
- Automatic internal incubation

Built-in Thermal Printer

- Real-time auto-printing
- Support for external printers

Designed Sampler

- Specially-designed sampler for small volume of specimen

Large Display

- 10" LCD Touch Screen
- Intuitive user interface

Connection

- Direct connection to LIS/HIS

Cartridge Tray


- Cartridges will be automatically ejected into the tray when test is finished

Test Reagent


- Room temperature storage
- Suitable for whole blood, serum, plasma and urine
- Get results in 3-15 minutes




TEST PROCEDURE




Sampling



Mixing



Dropping



Reading

APPLICATION



Laboratory



ER



ICU



Clinics



GP Office