

FinecareTM

**VITAMIN D RAPID
QUANTITATIVE TEST**

The era of high accuracy
in one-step testing

Wondfo

25-OH-VD

more than just a vitamin

New



One step - One step sample&buffer mixing



10 min - Get results in only 10 minutes



15 µL - Only 15 µL sampling volume




4-30°C - Room temperature storage




Whole Blood

Supports whole blood/plasma/serum sampling



Experience Innovation with Finecare™ Vitamin D Test

Simplified, faster, and more reliable, our brand new Vitamin D Reagent offers faster results with fewer steps and minimal sample requirements. Achieve accurate and reliable testing in just 10 minutes using whole blood—making VD testing easier than ever.



Vitamin D's Role in Health



Bone Health:

Rickets, Osteoporosis, Osteomalacia



Immune System:

Autoimmune, Infection Prevention



Maternal and Child Health:

Pre-eclampsia



Chronic Disease:


Diabetes, Cardiovascular

ITEMS SPECIFICATION


Item	Sample Type	Sample Volume	Loading Volume	Reaction Time	Detection Range	Storage
Vitamin D	Whole Blood/ Serum/Plasma	15 μ L	75 μ L/4 drops	10 minutes	5.0-100 ng/mL	4~30°C

TEST PROCEDURE


Regular Procedure




Sampling



Mixing




Loading



Testing

Capillary Sample Collector



Finecare Vitamin D is available on

1. Finecare™ FIA Meter Plus (FS-113)
2. Finecare™ II FIA Meter Plus SE (FS-114)
3. Finecare™ III FIA Meter Plus (FS-205)

