

FS-2000 Fluorescence Immunoassay Analyzer >



Specification

- Dimension: 325*475*395 mm
- Weight: 19 kg



Connectivity

- WIFI connection to Bi-directional LIS/HIS
- RS232 and ethernet port
- 5 USB ports



Built-in Code Scanner

- No need for ID chip
- Sample batch information



Built-in Waste Bin

- Automatic storage of used cartridges
- Holds up to 100 used cartridges



Special cooling function

- Ensures stability under high temperature



Multi-channel

- Semi-Automatic system
- 16 items tested simultaneously



Large Display

- 10-inch capacity touch screen
- Intuitive user interface
- Smooth operation experience like smart-phone



Built-in Thermal Printer

- Real-time auto printing
- External printer supported



Linux System

- Smooth operation
- Security system
- Custom intelligent management module
- Multi-user management

Operation Procedure >

1



Entry Test Information

2



Sample Loading

3



Start Testing

Technical Advantages >



Stability

Cutting edge TRFIA methodology
High stability, lower CV



Convenient

Simple operation procedure
Automatic collection of used cartridges



Instant Results

Results delivered in less than
8 minutes



High Efficiency

Maximum simultaneous detection
of up to 16 items



Delicate Design

Automatic reagent door
Charging indicator light
Forward USB port

1 QR Code Technology


- No need manual input cartridge information
- Reagent card recognition success rate is as high as 99%

2 Smaller Cartridge

- Lower shipping costs

3 Excellent CV and Accuracy with Optimized FIA Methodology

- Satisfying and stable results delivered in short time



Application Scope >



Ambulance



Emergency Room



Clinics



Pharmacy



GP Office